

ALPHA
TECHNOLOGIES


ARC 2020

Capillary Rheometer

**Direct measurement of
high shear rate viscosity**

With accurate die swell data

- Advanced material characterisation for Extrusion and Injection Moulding
- Invaluable for Process Development, Mould Design and Mould Flow




The First Name in Testing Solutions

Overview

The ARC 2020 Capillary Rheometer has been specifically designed for the Rubber and Thermoplastic Rubber Industries, and provides direct measurement of high shear rate viscosity, appropriate for QC/QA and development testing

- 12.7mm dia barrel for easier loading
- High piston force 15kN max
- Advanced material characterisation for Extrusion and Injection Moulding
- Ideal for QC, QA and Development
- Direct measurement of viscosity and die swell at high shear rates
- Shear Rate range: 1 – 29,000 s⁻¹
- Laser micrometer for die swell measurement
- Viscosity measurement over a 2 decade shear rate range in 10 minutes
- Invaluable for Mould Design, Mould Flow, and Process Development
- Temperature Range
50°C - 350°C (122°F - 662°F)



ARC 2020 benefits

The major shaping operations carried out in rubber processing include:

- Calendering
- Milling
- Extrusion
- Molding
 - Compression
 - Transfer
 - Injection

To optimize the processing conditions for a given batch of feedstock, the processability of the compound needs to be determined prior to processing.

Alpha Technologies have developed the ARC 2020 to provide rubber processors with the information required to optimize the processing conditions, for a specific production batch, thus maximize process efficiency and minimize scrap and process down time.

Major features

Flexible design


Interchangeable dies allows performance to be tailored to the application

- Die Diameter: 0.75mm to 2.0mm
- Die Length/Diameter Ratio: 15:1 to 30:1

Application	Die Diameter	Shear Rates*
Injection Molding	0.75 mm	2 to 29,000 1/s
Transfer Molding	1.00 mm	1 to 14,000 1/s
Extrusion	1.50 mm	0.2 to 4,000 1/s
Mixing	2.00 mm	0.1 to 1,700 1/s


Effect of Increasing Shear Rate

Increasing shear rate can differentiate between batches of compound. As can be seen below, two batches of compound show very different response to increased shear. The differing viscosity resulting from the increased shear will cause the compounds to 'process' differently.




Viscosity & Die Swell vs Shear Rate

Viscosity reduction with increased shear rate (shear thinning) influences the processability of compounds. The ARC2020 differentiates 'good' and 'bad' compound in a simple, quick test.



Die Swell Effects

Die swell can be measured under running and relaxed conditions. Running die swell influences the speed of extrusion/calendering process, while relaxed die swell predicts component dimensions.



Sample smoothness predicts component finish.



Good Fair Poor

LAB KARS for Windows™ Software

LAB KARS software is a powerful and easy to use rheological Windows™ based software package. It can easily accommodate routine calculations and plotting requirements for rheological analysis of capillary rheometer tests of rubber and thermoplastic materials. Alpha Technologies software can be used to quickly identify viscosity variations in batches and lots of materials. The program features Bagley and Rabinowitsch corrections, Carreau, Modified Cross, Power Law and Polynomial curve fits, Statistical error estimations, Shear rate dependence, Time at temperature relationships, critical shear stress, zero shear viscosity, and Intrinsic viscosity calculations.

LAB KARS Software for control and analysis

- Configure tests; Calculate results
- Plot rheological charts
- Bagley, Rabinovitch corrections
- Includes theoretical rheological models (eg. Carreau)

Specifications

Standards	ASTM D5099
Barrel	l = 229 mm (9.0") Dia. = 12.70 ± 0.005 mm (0.500 ± 0.0002")
Drive System	DC Servomotor
Piston Speed	0.03 to 600 mm/min., Resolution 0.01 mm/min. (0.001 to 23.6 in./min.)
Crosshead Force	15 kN Standard
Force Measurement	Load Cell, Piston; Barrel mounted transducer, automatic calibration
Piston Force Range	0 – 15 kN ± 0.5% (0 – 3370 lbs.)
Barrel Pressure Transducer	0 – 70,000 kPa ± 0.5% (0 – 10,000 psi)
Capillary Dies	Standard, (one) Stainless Steel die from list below: 0.75 mm (0.030"), 30:1 L/D 1.00 mm (0.040"), 30:1 L/D 1.50 mm (0.060"), 20:1 L/D 1.50 mm (0.060"), 15:1 L/D 2.00 mm (0.080"), 16:1 L/D
Die Swell Measurement	Laser Micrometer (± 0.03 mm)
Temperature Range	50°C - 350°C ± 0.2°C (122°F - 662°F)
Electrical	115 VAC ± 10%, 60 ± 3 Hz, 10 amp Single Phase 230 VAC ± 10%, 50 ± 3 Hz, 5 amp Single Phase
Data Processing System	LAB KARS for Windows™, includes PC with Windows XP
Dimensions	Width 47 cm (18.5 in.), Depth 56 cm (22 in.), Height 195 cm (76.8 in.)
Weight	Net 136 kg (300 lb) Gross 160 kg (350 lb)

Summary of features

- **Unique rubber testing features:**
 - Easy loading 12.7mm (0.5") dia. barrel
 - Force, pressure transducers at piston top, die entrance
 - High maximum crosshead force = 15 kN (3,370 lbs.)
 - High maximum barrel pressure = 70,000 kPa (10,000 psi)
 - LAB KARS software for analyzing both rubber and thermoplastic test results
 - Extra length 229 mm (9") barrel for preheating rubber
- **Capable of measuring viscosity over a 2 decade shear rate range in a multi-zone test of less than 10 minutes:**
 - 100 – 10,000 s⁻¹ with a 0.75 mm (0.030") dia. die
 - 30 – 3,000 s⁻¹ with a 1.0 mm (0.040") dia. die
 - 10 – 1,000 s⁻¹ with a 1.5 mm (0.060") dia. die
 - 3 – 300 s⁻¹ with a 2.0 mm (0.080") dia. die
 - Wider shear rate range available with longer test times

For additional information on Alpha Technologies products please contact:

Alpha Technologies Services LLC

3030 Gilchrist Road
Akron, Ohio 44305, USA
Tel +1 800 356 9886
Fax +1 330 848 7326

Alpha Technologies UK

Unit 2B Crowood House
Gipsy Lane, Swindon SN2 8YY, UK
Tel +44 1793 601 100
Fax +44 1793 615 214



ALPHA
TECHNOLOGIES

 **Dynisco**
company

Disclaimer of warranty and liability

The information provided herein is believed to be true and correct but no warranty is given to its completeness, accuracy or fitness for use for any particular purpose. Alpha Technologies shall not be liable for any loss or damage arising from the failure to achieve particular result by the application of any route, method or process that is recommended herein.

www.alpha-technologies.com