

# **COMPANY PROFILE**

## **CONTENTS**



ABOUT US

**VISION** 

PRODUCT OVERVIEW

**APPLICATIONS** 

SALES STRATEGY

**PRODUCTS** 

CLIENTS LIST

INTERNATIONAL PARTNERSHIPS

CONTACT US

## **ABOUT US**



### EQUALITY has two meanings, Excellent quality and Equality among all employees.

EQUALITY Co. aims to provide equipment to the customer pursuing the best quality, to maintain the principle of equality among the employees to live well together.

We strive to become a company specializing in whole equipment for research and experimentation in ceramic, chemical, pharmaceutical and materials fields and pilot plant facilities.

We supports development and testing application technologies for equipment used in existing university laboratories and corporate laboratories.

And through this, we will surely supply the best equipment for various experiments to the customer!



## **VISION**





### **PRODUCTS OVERVIEW**



EQUALITY Co. is dealing with whole laboratory and industrial equipment for-MIXING, DISPERSING and MILLING samples.

### PASTE TYPE



Paste syringe charger

### **DRY TYPE**



### WET TYPE



### **ANALYSIS & OTHERS**



\* Please refer to our general catalog for specific information

### **APPLICATION**



**EQUALITY Co.** contacts customers who in variety of industrial fields such as,





Automobile





**Battery** 

















As we are dealing with basic and essential equipment in material processing, We keep exploring various clients in tremendous industry!

## **SALES STRATEGY**



### Counselling

Visiting customer side for meeting

Recommending suitable equipment for the application.

Operating demo test room (Whole demo machines equipped)



### Wide-Sales channel

Participating exhibition (7 to 10 times/a year)

**Direct contact** 

Web-site, Blog, Youtube

### **Technical capability**

Maintenance & repairing service (Skilled & experienced engineers)

Product customizing (for domestic products)

### Reserch

Exploring new items (matching customer's requirements)

Visiting exhibition (Local/Int'l)

Documentation of test applications










#### Three Roll Mill

Equipment for dispersing paste samples using shear force between rollers generated by rotating the three rollers in different directions.



- ✓ Control 0.1mm to adjust 0.1µm gap size
- √ High resolution gap adjustment
- ✓ Easy visible gap size displayed







LCD display

Various roller material

<sup>\*</sup> Please refer to our general catalog for specific information



#### PASTE MIXER / INTER MIXER

Paste mixer/Inter mixer is using planetary rotation and revolution to mix and defoam the sample cleanly and quickly without impeller. High capacity equipment capable to produce up to 20kg. Dedicated adapter for customer's own container and vial



**PDM-300V** 



IMX-500CV




PDM-20KV





Centrifugal rotation of paste mixer/inter mixer.









Cooling sample function for Inter mixer





Wireless temp. sensor/Display



Adapter customizing

<sup>\*</sup> Please refer to our general catalog for specific information



#### SYRINGE CHARGER

For charging paste samples into packaging syringe. Possible to charge both low and high viscosity sample.

Throttle valve type and Screw extruding type



Throttle method

For low viscosity samples(Throttle type)





For high viscosity samples(Screw type)











1Step: Pressing

2Step: Screw Extruding

3Step: Charging









Dental Material Paste Charging (Viscosity: About 1,000,000cps)



#### POWDER MIXER / Air Dry Mill / DISK MILL / JET MILL

Dry type processing equipment for mixing and milling of powdery samples.



Air Dry Mill(ADM-5)











Powder mixer









Air Jet Mill



Disk mill









<sup>\*</sup> Please refer to our general catalog for specific information



#### KM MILL / SUPER MILL / ATTRITION MILL / BASKET MILL

Wet type processing equipment for mixing and dispersing of slurry, liquid samples.












**Basket mill** 



<sup>\*</sup> Please refer to our general catalog for specific information



#### Particle size analyzer / Table-top type SEM

We offer Particle size analyzer and SEM for analyzing equipment after the production as below.









**CUBE II** 

<sup>\*</sup> Please refer to our general catalog for specific information



#### Rotary evaporator / Ultrasonic disperser / Basic lab equipment

You can meet other various equipment for mixing/dispersion/milling etc.



<sup>\*</sup> Please refer to our general catalog for specific information

### **CLIENTS LIST**



### **Electronics / Materials**













































Semisysco



Sammi

















INTEROJO Lict







### Pharmaceuticals & Cosmetics















CWM(주)코웰메디 🛟 NeoPharm





















### Universities & Ohters



**)** 영남대학교
























### **INERNATIONAL PARTNERS**



EQUALITY Co. is engaged in the import/export of industrial equipment and sales marketing based on its excellent technology and overseas sales capability.



## **CONTACT US**



#### **ADDRESSES**

- Office/Lab: #403 Creation-hall, 61 Dongbu-ro, Chein-gu, Yongin-si, Gyeonggi-do, Republic of Korea(17145)
- Factory: 420-9, Jungbu-daero, Majang-myeon, Icheon-si, Gyeonggi-do, Republic of Korea(17385)


### TEL / FAX

- Tel: +82 031 333 2346(Office)

- Fax: +82 031 335 2346

#### WEBSITES

- www.ex-quality.com
- E-mail: e-quality@naver.com, equalitywon@naver.com
- Youtube: https://www.youtube.com/channel/UCnhkNO0yctd-c-pDIZq9T7Q





# **THANK YOU!**