Generation of Ultrashort Pulses

introducing Series

NEW!

Ultrafast Harmonic Generator for Femtosecond and Picosecond Oscillators

GWU-Lasertechnik

Specifications

- high conversion efficiency for maximum harmonic power

Femtosecond	2ω (SHG)	3ω (THG)	4ω (FHG)
Wavelength [nm]	340 – 540 ⁽¹⁾	226 – 360 ⁽¹⁾	209 – 242 (2)
Conversion Efficiency	40%	10%	4%

Picosecond	2ω (SHG)	3ω (THG)	4ω (FHG)
Wavelength [nm]	340 – 540 ⁽¹⁾	226 – 360 ⁽¹⁾	209 – 242 ⁽²⁾
Conversion Efficiency	20%	10%	on demand




⁽¹⁾ Input wavelength 680 nm – 1080 nm

⁽²⁾ Input wavelength 836 nm - 986 nm

Tuning characteristics

- ultrafast from the visible to the deep-UV



Typical tuning curves for femto operation with Spectra Physics Tsunami XP lasers


Input: 680 – 1080 nm

- 2ω (SHG):340 540 nm
- 3ω (THG)226 360 nm
- 4ω (FHG)
 209 242 nm



Ultrashort pulses

- femto and pico operation




Typical autocorrelation trace for femto operation, pump 80 fs

- minimum pulse broadening
- easy conversion between femto and pico operation



Excellent beam parameters




Typical beam profiles for femto operation



Reliable operation

- excellent long term stability



Typical power stability with Spectra Physics Mai Tai HP

