


**HANAROSTENT®**

**Duodenum/Pylorus Lasso (NCN)**



**Indication :**

- Palliative treatment of pyloric or duodenal stricture caused by malignant tumors

**Features :**

- Flare heads on the both ends for Anti-Migration
- Uncovered on the both ends for Anti-Migration
- Lasso on the proximal

**Endoscopic Application**


	Stent (mm)			Delivery device	
	Diameter	Usable length	Total length*	Length(mm)	Diameter(mm/fr)
NDCL-18-060-230	23-18-23	24	60	2300	3.4/10.2
NDCL-18-150-230		114	150	2300	3.4/10.2
NDCL-20-060-230	26-20-26	24	60	2300	3.4/10.2
NDCL-20-150-230		114	150	2300	3.4/10.2

\*Increment of 10mm

	(NNN) Non covered	<u>XXX-00-000-000</u>	<u>XXX 00 000-X000</u>
	(CCC) Fully covered	① ② ③ ④	① ② ③ ④
	(NCN) Partially covered	① Product Name	③ Stent Length
		② Stent Diameter	④ Delivery Device Length

**HANAROSTENT®**

**Duodenum/Pylorus Lasso (NCN)**



**Indication :**

- Palliative treatment of pyloric or duodenal stricture caused by malignant tumors

**Features :**

- Flare heads on the both ends for Anti-Migration
- Uncovered on the both ends for Anti-Migration
- Lasso on the proximal

**Application a fluoroscopic Application**

	Stent (mm)			Delivery device	
	Diameter	Usable length	Total length*	Length(mm)	Diameter(mm/fr)
NDCL-18-060-180	23-18-23	24	60	1800	3.4/10.2
NDCL-18-150-180		114	150	1800	3.4/10.2
NDCL-20-060-180	26-20-26	24	60	1800	3.4/10.2
NDCL-20-150-180		114	150	1800	3.4/10.2

\*Increment of 10mm

	Stent (mm)			Delivery device	
	Diameter	Usable length	Total length*	Length(mm)	Diameter(mm/fr)
NDCL-18-060-120	23-18-23	24	60	1200	3.4/10.2
NDCL-18-150-120		114	150	1200	3.4/10.2
NDCL-20-060-120	26-20-26	24	60	1200	3.4/10.2
NDCL-20-150-120		114	150	1200	3.4/10.2

\*Increment of 10mm

	(NNN) Non covered	<u>XXX-00-000-000</u>	<u>XXX 00 000-X000</u>
	(CCC) Fully covered	① ② ③ ④	① ② ③ ④
	(NCN) Partially covered	① Product Name	③ Stent Length
		② Stent Diameter	④ Delivery Device Length