AnyNeb SUPER MESH NEBULIZER

- Vibrating mesh technology
- No metal elution (Nickel & SUS Free)
- Noiseless operation



 AnyNeb is ideal for the On-The go treatment of lower airways respiratory diseases like asthma even in public places.

- AnyNeb sprays liquid medication into lung as aerosols.
- AnyNeb allows using various type of liquid medication directly without dilution.
- Vibrating mesh releases medication via 4um of fine holes.
- Noiseless operation.
- Residual volume is under 1% by world patented technology.
- Inhalation is possible at over 75° of tilt angle.





⊂€ 0120




model **NE-SM2**

AnyNeb SUPER MESH NEBULIZER



Technical Data



Type Nebulizer Type	Pulmicort	Ventolin
Ultrasonic Nebulizer	Bad	Available
Mesh Nebulizer	Excellent	Excellent

Color Option



Specifications

Power & Consumption Battery	AA Alkaline Battery 3V DC 2ea	
Dimensions	(W) x (D) x (H) = 47*65*106(mm)	
Weight Main Unit	88g, Mask : 50	
Nebulization Rate Approximately	0.25ml/min(23°C)	
Particle Size Average	5μm	
Volume of Medication Bottles	6ml	
Operating Time	>1 hour 30 minutes	
Residual Volume	<0.1ml	
Ultrasound Frequency Approximately	108KHz	
Operating Environment Temperature	10~40°C Humidity : 30~85% R.H	
Pressure	70~106KPa	
Storage Environment Temperature	-20~60°C Humidity: 10~95% R.H	
Main Unit	5 years from purchase	
Mesh Cartridge	1 year from purchase(condition: 3times/1day	
Instrument Warranty	1 year from purchase	
Contents	Nebulizer Unit, Adult Mask, Child Mask, Mouthpiece, Alkaline Battery AAx2(1.5V),	

*All specifications are subject to change prior notice.

Components

Mesh Cartridge	Medication Bottle	Mask Adapter
Mouthpiece	Small Inhalation Mask	Large Inhalation Mask

*For repairing or repurchasing of parts, please call our company or local distributor.



KTMED INC. #602, Byucksan Digital Vally 1-cha, 24, Digital-ro 27-gil, Guro-gu, Seoul, Korea **Tel**: +82-2-2108-7890 **Fax**: +82-2-2108-7894

